



## TLC-FlashReader™

Your multi-function TLC plate reader

Accepted TLC plate size	Any TLC plate size up to 8.5cm wide and 12.5cm high
Light sources	Illumination with white light (2 lamps), UV 254nm (2 lamps) and UV 366nm (2 lamps)
Use	"stand alone" mode: controlled by "TLC-FlashReader™" software allowing the visualization and archiving of TLC plates Or coupled with InterSoft® X software for direct transfer of the TLC plate into Genius and archiving of TLC plates
Safety	Automatic shut off for the lamps when the drawer is opened to protect the user from UV radiation
Cleaning	Removable TLC plate holder for easy maintenance
Compatibility	The TLC-FlashReader™ can be paired with all of our puriFlash® Gen 5, puriFlash® XS 520 Plus & puriFlash® 535-XL devices
Communication	2 USB cables
Power supply	Power cord of choice: China, India, UK, US, European & Switzerland Voltage: 220V or 110V
Dimensions & Weight	Weight: 5.9kg

