

Evaporator puriFlash® XS-Vap

**Select the right configuration (accessories + instrument)
of puriFlash® XS-Vap according to customer needs**



PART I : SELECT THE VOLTAGE OF THE PURIFLASH XS-VAP

Step 1 Select the voltage of your unit	Voltage 220-240V Delivered with (1) Basket ; (1) Gas inlet tubing 6mm (length : 1x3 meters); (1) Drain bath tube (length : 1x2 meters) P/N: EVAA00 (1u)			or	Voltage 100-120V Delivered with (1) Basket ; (1) Gas inlet tubing 6mm (length : 1x3 meters); (1) Drain bath tube (length : 1x2 meters) P/N : EVAA01 (1u)	
Step 2 Select the power cord according to the country	Europe P/N power cord : XZ2190 (1u)	UK P/N power cord : XZ2270 (1u)	Switzerland P/N power cord : EVAA90 (1u)	China P/N power cord : EVAB00 (1u)	India P/N power cord : EVAB10 (1u)	US P/N power cord : XZ2340 (1u)

PART II : SELECT THE ACCESSORIES OF THE PURIFLASH XS-VAP ACCORDING TO THE CUSTOMER REQUEST

For small sample volumes (up to 25mL) - Test tubes and vials					or / and	For high sample volumes (up to 250mL) - glassware 250mL capacity
Step 3 Select the manifold head	Manifold 56 positions for test tubes 12,13,16,18mm OD (like 2mL vials, 13x100 & 16x100mm test tubes) Manifold with 4 rows of 14 positions delivered with the grid Up to 56 samples in parallel Delivered without nozzles (p/n : EVAA50) P/N Manifold : EVAA20 (1u)			Manifold 90 positions for test tubes 13mm OD (like 13x100mm test tubes) Manifold with 5 rows of 18 positions delivered with the grid Up to 90 samples in parallel Only suitable with rack P/N : 1R8590 (1u) Delivered without nozzles (p/n : EVAA30) P/N Manifold : EVAA30 (1u)		Manifold for 250mL glasswares Up to 6 samples P/N Manifold : EVAA40 (1u)
Step 4 Select the nozzle	Nozzles for manifold EVAA20 & EVAA30 If you use manifold EVAA20, you can offer up to 4x EVAA50 (either 4x14 nozzles) If you use manifold EVAA30, you can offer up to 7x EVAA60 (either 7x14 nozzles) Nozzles can be disconnected manually from manifold. You can also use the tool p/n : EVAA60 P/N Nozzles : EVAA50 (pack of 1x14u)					No nozzles required Nozzles (6) (p/n : 390191) are already included with manifold EVAA40
Step 5 Select the adapter	Adapter for rack test tubes Only if you use manifold EVAA20 and EVAA30 The adapter must be connected inside the basket delivered with puriFlash XS-Vap If the customer wants to evaporate 18x150mm tubes with a temperature >60°C, please offer the option p/n : EVAA80 P/N adapter : EVAA10 (1u)					No adapter required No adapter needed with the manifold 6 positions (p/n : EVAA40)
Step 6 Select the rack according to the customer requests	Rack especially designed for GC / LC users All the racks are positioned in the adapter EVAA10 The manifold suitable with these racks is EVAA20 Not suitable with the flash chromatography collector Interchim					
	RACK 56 POSITIONS_15ML_16X100MM TEST TUBES_FOR PURIFLASH XS-VAP P/N rack : EVAB30 (1u)	POSITIONS_10ML_13X100MM TEST TUBES_FOR PURIFLASH XS-VAP P/N rack : EVAB40 (1u)	POSITIONS_2ML_12X32MM VIALS_FOR PURIFLASH XS-VAP P/N rack : EVAB50 (1u)			
	RACK 56 POSITIONS 8ML 17X60MM VIALS FOR PURIFLASH XS-VAP P/N rack : EVAB80 (1u)	RACK 56 POSITIONS_2ML_CENTRIFUGE TUBES_FOR PURIFLASH XS-VAP P/N rack : EVAB60 (1u)	RACK 56 POSITIONS_4ML_15X45MM VIALS_FOR PURIFLASH XS-VAP P/N rack : EVAB70 (1u)			
	Racks used with the puriFlash chromatography systems Interchim (Gen4, Gen5 and XS) All the racks are positioned in the adapter EVAA10 The manifold suitable with these racks is EVAA20			Racks used with the puriFlash chromatography systems Interchim XS The rack is positioned in the adapter EVAA10 The manifold suitable with this rack is EVAA30		
	16X150 MM RACK Rack for PF215 , PF430, PF450 , PF4125, PF425 P/N rack : PF4381 (1u)	18X150 MM RACK rack for PF-XS520, PF-XS535 P/N rack : 1R8610 (1u)	18X150 MM RACK Rack for PF215 , PF430, PF450 , PF4125, PF4250 P/N rack : PF4391 (1u)	13X100 MM RACK rack for PF-XS520, PF-XS535 only suitable with manifold p/n : EVAA30 P/N rack : 1R8590 (1u)		
	16X150 MM RACK Rack for PF5020, PF5050, PF5125, PF5250, PF5250P P/N rack : AYHE40 (1u)	18X150 MM RACK Rack for PF5020, PF5050, PF5125, PF5250, PF5250P P/N rack : AYEH50 (1u)				
Step 7 Select the glasswares	Please check the list of our test tubes					GLASSTUBE 250ML WITH 1ML END-POINT_OD 56MM_THICKNESS 2.5MM P/N glassware : BVLGH0 (1u) CENTRIFUGE TUBES_56,00X147MM_BOROSILICATE GLASS DURAN_250ML_ROUND BOTTOM P/N glassware : B684D0 (1x10u)
Step 8 Select the fume enclosure if the customer wants to use the system outside of fume hood	FUME ENCLOSURE WITH VENTING HOSE 4 METERS The fume enclosure is connected on the rear of the puriFlash XS-Vap Suitable for puriFlash XS-Vap EVAA00 (220-240V) and EVAA01 (100-120V) P/N fume enclosure : EVAA70 (1u)					

[illegible][illegible][illegible][illegible][illegible][illegible]